

# HIGH PERFORMANCE LIQUID CHROMATOGRAPHY (HPCL)


*Paper III CA 2: Interpretation of HPCL Chromatogram*

<b>Name</b>	Hansa Sharma
<b>Roll No.</b>	52
<b>UID</b>	22BSC159

# Interpretation of given HPCL Chromatograms

**Standard: Aspartame**

**Unknown: Coke Zero in Aspartame**



Aspartame Samples	Peak Height (Absorbance)	Peak width (seconds)	Concentration (mg/L)
Standard 1	0.267	12.11	122.2
Standard 2	0.487	11.76	244.4
Standard 3	0.847	11.25	366.6
Standard 4	1.070	11.85	488.8
Standard 5	1.250	12.98	611.0
Unknown (Coke Zero)	0.386	11.25	-

## **Peak analysis**

1. The standard chromatogram displays clear peaks for each aspartame concentration.
2. The aspartame peak in Coke Zero is distinctly separated.
3. There is no co-elution detected in either chromatogram.
4. The standard chromatogram shows background noise between 50 to 100 seconds.
5. The chromatogram for the diet soda exhibits significant noise from approximately 40 to 100 seconds.
6. The aspartame peaks appear as skewed, asymmetric curves.
7. Both chromatograms maintain stable baseline readings.
8. The peaks in the standard chromatogram are well-defined for each concentration of aspartame.

## **Approximate Retention time**

1. Aspartame Standard 1= 185 seconds
2. Aspartame Standard 2= 183 seconds
3. Aspartame Standard 3= 180 seconds
4. Aspartame Standard 4= 182 seconds
5. Aspartame Standard 5= 184 seconds
6. Aspartame in Coke Zero= 185 seconds

# Calculation of Peak Area

## Peak area for Standard 1

- Peak Height (PH): 0.267 (mm)
- Peak width at half height (PW): 12.11s
- Peak Area =  $PH \times PW$   
 $= 0.267 \times 12.11 = \mathbf{3.23 \text{ mm}^2}$

## Peak area for Standard 2

- Peak Height (PH): 0.487 (mm)
- Peak width at half height (PW): 11.76s
- Peak Area =  $PH \times PW$   
 $= 0.487 \times 11.76 = \mathbf{5.73 \text{ mm}^2}$

## Peak area for Standard 3

- Peak Height (PH): 0.847 (mm)
- Peak width at half height (PW): 11.25s
- Peak Area =  $PH \times PW$   
 $= 0.847 \times 11.25 = \mathbf{9.53 \text{ mm}^2}$

## Peak area for Standard 4

- Peak Height (PH): 1.070 (mm)
- Peak width at half height (PW): 11.85s
- Peak Area =  $PH \times PW$   
 $= 1.070 \times 11.85 = \mathbf{12.69 \text{ mm}^2}$

## Peak area for Standard 5

- Peak Height (PH): 1.250 (mm)
- Peak width at half height (PW): 12.98s
- Peak Area =  $PH \times PW$   
 $= 1.250 \times 12.98 = \mathbf{16.23 \text{ mm}^2}$

## Peak area for Coke Zero

- Peak Height (PH): 0.386 (mm)
- Peak width at half height (PW): 11.25s
- Peak Area =  $PH \times PW$   
 $= 0.386 \times 11.25 = \mathbf{4.34 \text{ mm}^2}$

## Calculation of Percentage

Sample	Peak Area of Standard 1	Peak Area of Standard 2	Peak Area of Standard 3	Peak Area of Standard 4	Peak Area of Standard 5	Peak Area of Coke Zero
<b>Caffeine</b>	2.2	4.0	7.27	9.56	12.08	9.04
<b>Aspartame</b>	3.23	5.73	9.53	12.69	16.23	4.34
<b>Benzoate</b>	3.44	5.49	9.88	14.5	18.06	15.18

### Percent Aspartame in Coke Zero

$= (\text{Area of Aspartame curve for Coke Zero} \div \text{Area of the total curve}) \times 100$


$= (4.34 \div (3.23 + 5.73 + 9.53 + 12.69 + 16.23)) \times 100$

$= (4.34 \div 47.41) \times 100$

$= 9.154 \%$

# Estimation of Aspartame Concentration in Coke Zero

The calibration curve for Aspartame



Peak Area of Aspartame *for* Coke Zero = 4.34mm<sup>2</sup>

Peak Area of Aspartame Standard 1 = 3.23mm<sup>2</sup> (This is the closest to the unknown peak)

The Concentration of Aspartame in Standard 1 is 122.2 mg/L

∴ The concentration of aspartame present in as is equal to 122.2 mg/L

But the Diet Soda was diluted by a dilution factor of 2. Thus, the actual concentration of Aspartame in would be:

$$122.2 \text{ mg/L} \times 2 = 244.4 \text{ mg/L}$$